

Innovations for Quality in Life

The Smart Solution for all Medical Bags

SmartBagEco uses flexible sensors to measure the fill level of medical fluid bags. As the first solution developed for this application, it combines an advanced measurement principle with user-friendly design, making it accessible to everyone. This ease of use supports nursing tasks and enhances patient care. The system is compatible with all types of fluid bags, such as stoma, infusion and urine bags, making it versatile across applications in hospitals, outpatient care, and even private use at home.



SmartBag
Ecosystems


- Compatible with **all medical fluid bags**. ✓
- Scalable fill level monitoring. ✓
- Personal and professional use. ✓
- Continuous fill level measurement. ✓

Our Ecosystem

Our Ecosystem essentially consists of three components: *the Sensor, the Clip-On Module* and *App*. Designed for usage with all medical Bags. **Made for Carers, Patients** and **Family Members**.



The Sensor measures the fill level and simply needs to be attached to the medical fluid bag via an adhesive patch. Our sensor is fabricated by the use of **screen printing**, making it a time and cost efficient way of produce. It is also integrated within an adhesive bond patch for easy handling. It just needs to be attached to the medical bag currently in use.
hygiene focused single use item



The Clip-On Module features a 3D-printed casing that houses the electronics and incorporates a spring mechanism for secure and stable contact with the sensor. This design ensures reliable data transfer and includes alignment points for easy handling and accurate placement, simplifying setup for the user.
reusable item



The App offers customizable settings, allowing users to *adjust bag type, volume, and fill level thresholds*. Its intuitive design features large icons and simple navigation for ease of use. Built with **Flutter**, the app is *cross-platform ready*, ensuring adaptability on iOS, Android, and Windows. Utilizing **Bluetooth Low Energy (BLE)**, it provides reliable, battery-friendly communication with the sensor.
digital companion



Value for all Stakeholders




Measurement with effusivity

We measure the fill level based on **effusivity (e)**, which is defined as:

$$e = \sqrt{k \cdot \rho \cdot c}$$

using the relation $\Delta\vartheta(t) = \frac{2q}{e\sqrt{\pi}} \sqrt{t}$



Learn more about our SmartBagEco. Scan the QR Code below or use our NFC tag.

